

Plasticity

Volumetric Specimen Cutter

Part NO: WAS2

- **Easy to use**
- **Higher productivity**
- **Compact robust design**

This new Specimen Cutter with its compact robust design allows the user to easily cut multiple volumetric samples for use with the P14 Plastimeter and O14 Ageing oven. An easy to operate circular motion of the handle, allows the cutting of consistent volumetric samples quickly and easily.



Volumetric Specimen Cutter	
Weight	3.8 kg
Dimensions (excl handle)	W 10cm D 11.5cm H 17cm
Standards	ISO 2007 ASTM D3194

